LABORATÓRIÁ CENTRA STU PRE NANODIAGNOSTIKU

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
LABORATORIES OF STU CENTRE FOR NANODIAGNOSTICS

Description of main activities:

Analysis of micro/nanostructures and systems for micro/nanoelectronic applications, sensorics, photonics, electronics, medicine, geology, material engineering and environment. Analysis of interfaces is possible at subangstrom level. Devices provide a possibility to examine:

- chemical composition,
- crystallography, defects, and crystallographic orientation,
- manifestations of properties of electron structure,
- surface contamination,
- 3D morphology of objects in nanoscale,
- microscopy of magnetic materials,
- optimisation of technological processes.

At present, the Centre deals with research of interaction of volume and 2D carbon materials, as they are promising and prospective materials for electronic applications.

Equipment of the Transmission Electron Microscopy Laboratory

Analytical transmission electron microscope JEM ARM 200cF with the resolution of 0.78 Å equipped with field auto emission electron source and with

two correctors of spherical aberration enables a complex characterization of nanostructured objects on the subatomic level by means of the following methods: imaging with high angle annular dark field detector (HAADF), imaging with low angle annular dark field detector (LAADF), imaging with annular bright field detector (ABF) in scanning regime (STEM), in TEM/HRTEM regimes, and imaging in back-scattered electron (BEI) and secondary electron regimes (SEI). Electron diffraction methods (selected area diffraction, nanobeam diffraction, convergent and precession diffraction) are available to acquire the information related to crystallographic structure. Further possibilities include orientation and phase mapping, qualitative and quantitative energy dispersive analysis (EDS), electron energy loss spectroscopy (EELS) and mapping of atom distribution in the sample. Microscope has incorporated Lorentz mode for characterization of magnetic samples and biprism for electron holography. 3D morphology of nanoscaled objects can be obtained by employing electron tomography. Other adjacent laboratories serve for preparation of TEM samples where various physical and chemical methods are exploited.

Equipment of the Auger Spectrometry Laboratory



The Auger microprobe JEOL JAMP 9510-F serves for elemental analysis of various structures and it combines the advantages of scanning electron microscopy and Auger electron spectroscopy. Except of standard Auger electron spectra and concentration depth profiles, recording of Auger maps and line profiles is possible. In secondary electron mode, the best resolution is 3 nm and the smallest spot of the electron beam in the Auger analysis mode is 8 nm. Auger spectroscopy is an appropriate complement of methods incorporated in TEM, because Auger signal is released from the topmost surface of a structure under study, usually from 0.5 – 3 nm depths. In the mode with energy resolution $\Delta E/E = 0.05\%$ exists a possibility to examine also plasmonic excitations. Owing to these parameters, the device is dedicated to analysis of nanostructures.

CONTACT

Slovak University of Technology in Bratislava


Vazovova 5
812 43 Bratislava
Slovak Republic

Assoc. Prof. Ing. Marian Veselý, CSc.
marian.vesely@stuba.sk

STEM image of graphene monolayer and multilayer and EDS map of Ga atoms (green) and As (red) in GaAs

A part of transmission electron microscope Jeol JEM ARM 200cF

Analytical chamber of Auger microprobe Jeol JAMP 9510-F